

Respiratory Fit Testing in Agricultural Communities

...HELPING FIGHT THE ENEMY WE CANNOT SEE

STATE OF SCIENCE - 2020



SEPTEMBER 18, 2020

AGRISAFE NETWORK

CHARLOTTE HALVERSON BSN,COHN-S


CLINICAL DIRECTOR, AGRISAFE NETWORK



CDC/NIOSH partnership




NIOSH Centers for Agricultural Safety and Health



- | | |
|---------------------------------------|--|
| 1. Bassett Healthcare | 6. University of Iowa |
| 2. Colorado State University | 7. University of Kentucky |
| 3. National Farm Medicine Center | 8. University of Minnesota |
| 4. University of California, Davis | 9. University of Nebraska Medical Center |
| 5. University of Florida, Gainesville | 10. University of Texas Health Science Center, Tyler |
| | 11. University of Washington |

Safety Considerations – What Can We Do?

Engineer dangers out of environment if possible!



Personal protective equipment - may be only feasible solution
The last line of defense!


Know Your Exposures – There are Many!



Gas	Health Effects		
	Low	Medium	High
Hydrogen Sulfide (H₂S)	2-20 ppm: nausea, headache, dizziness	100-300 ppm: altered breathing, fluid in lung	500-700 ppm: collapse, death
Methane (CH₄, LEL)	5-15% (100% LEL): asphyxiant	< 1% (10,000 ppm): No known toxicity	< 5% (50,000 ppm): Explosive
Ammonia (NH₃)	5-20 ppm: odor, eye irritation	20-50 ppm: Moderate eye and upper respiratory tract irritation	2500 ppm: chemical pneumonitis, edema, cyanosis, death
Carbon Dioxide (CO₂)	600-1000 ppm: odor, muscle stiffness	1000-2500 ppm: general drowsiness; 5000 ppm: 8-hr maximum	30,000 ppm (3%): increased pulse rate, nausea, impairment
Carbon Monoxide (CO)	< 9 ppm: comfortable living concentration (35 ppm = 8-hr allowable)	200 ppm: headache, dizziness, nausea in 35 ppm = 8-hr 2 hours	400 ppm: life threatening in 3 hours

Children, elderly, pregnant women, and other risk factors: CO, H₂S, and NH₃. The concentrations are advisory only as "no-harm" levels.





When wearing a respirator – the fit is the key

- ▶ All respirators - N95 and above – must fit tight to the face and allow no air leaks
- ▶ If worn correctly, a respirator will filter out 95% or more of aerosol/particulates
- ▶ Respirators must:
 - ▶ be donned correctly
 - ▶ have a fit check performed each time it is worn
 - ▶ be removed (doffed) properly






Image courtesy of the U.S. EPA

Particle Size

Seal Check (Fit check) Procedure: an employee (wearer) responsibility – and is done each time a respirator is used

POSITIVE PRESSURE CHECK



NEGATIVE PRESSURE CHECK



Very Important!

Seven Steps to Correctly Wear a Respirator at Work

Following these simple steps will help you properly put on and take off your respirator, and keep you and everyone else safe.

1 Wash Your Hands




Wash your hands with soap and water or alcohol-based hand rubs containing at least 60% alcohol.

2 Inspect the Respirator



Inspect the respirator for damage. If it appears damaged or damp, do not use it.

3 Put on the Respirator



Clip the respirator in your hand with the nosepiece at your fingertips and the straps hanging below your head.



Cover your mouth and nose with the respirator and make sure there are no gaps (e.g., facial hair, hair, and glasses) between your face and the respirator.



Place the strap over your head and rest at the top back of your head. If you have a second strap, place the bottom strap around your neck and below your ears. Do not crisscross straps.



If your respirator has a metal nose clip, use your fingertips from both hands to mold the nose area to the shape of your nose.

4 Adjust the Respirator



Place both hands over the respirator. Inhale quickly and then exhale. If you feel leakage from the nose, readjust the nosepiece; if leakage from the respirator edges, readjust the straps.



Repeat until you get a proper seal. If you can't get a proper seal, try another respirator.



Avoid touching the respirator while using it. If you do, wash your hands.

Note: If you reuse your respirator, wear gloves when inspecting and putting on the respirator. Avoid touching your face (including your eyes, nose, and mouth) during the process.

6 Remove the Respirator



Wash your hands.



Remove the respirator from behind. Do not touch the front.



If the respirator does not need to be reused because of supply shortages, discard it in a closed-bin waste receptacle. Wash your hands.

For more information, see the quick video, "Putting On and Taking Off a Mask".



Helping You Wear it Right

Wearing Your Filtering Facepiece Respirator

1



Place the respirator over your nose and mouth. Be sure the metal nose clip is on top. With models 8210 or 07048, pre-stretch the straps before wearing.

2



Pull the top strap over your head until it rests on the crown of your head above your ears.

3



Pull the bottom strap over your head until it rests just below your ears.

4



Using both hands starting at the top, mold the metal nose clip around your nose to achieve a secure seal.



Filtering Facepiece (non-valved) Respirator

Filtering Facepiece (non-valved) Respirator

Check the Seal of Your Filtering Facepiece Respirator Each Time You Don the Respirator.



Positive Pressure User Seal Check

For Non-Valved Respirators

Place both hands completely over the respirator and *exhale*. The respirator should bulge slightly. If air leaks between the face and faceseal of the respirator, reposition it and readjust the nose clip for a more secure seal. If you cannot achieve a proper seal, **do not** enter the contaminated area. See your supervisor.



Negative Pressure User Seal Check

For Valved Respirators

Place both hands over the respirator and *inhale* sharply. The respirator should collapse slightly. If air leaks between the face and faceseal of the respirator, reposition it and readjust the nose clip for a more secure seal. If you cannot achieve a proper seal, **do not** enter the contaminated area. See your supervisor.

WARNING

Respiratory Products

These respirators help protect against certain airborne contaminants. Before use, the wearer must read and understand the User Instructions provided as a part of the product packaging. A written respiratory protection program must be implemented meeting all the requirements of OSHA 1910.134 including training, fit testing and medical evaluation. In Canada, CSA Standard Z94.4 requirements must be met. Respiratory requirements of the applicable jurisdiction, as appropriate. **Failure may result in sickness or death.** For proper use, see packaging instructions, supervisor, or call 3M Chemical Technical Service in USA at 1-800-243-6030 and in Canada at 1-800-387-4434.

General Office:
3M Center
St. Paul, MN 55144-1000

Please recycle. Printed in U.S.A.
© 2014 3M. All rights reserved.
70-0707-0700-0



OSHA Occupational
Safety and Health
Administration
www.osha.gov

1-800-321-OSHA (6742)
TTY 1-877-889-5627

OSHA 4015-02 2020

<https://www.osha.gov/Publications/OSHA4015.pdf>

<https://www.youtube.com/watch?v=oU4stQgCtV8>

Filtering Facepiece Fit Test

Fit Testing is required if:

- ▶ your employer requires the use of respirators.
- ▶ you (the employee) are exposed above the OSHA permissible exposure levels (PELs) of particulate or aerosol.



*If workers wear these respirators
on a voluntary basis NO fit testing
is required.*



Type of Fit Testing

- **Quantitative**

- Portacount test

- **Qualitative (pass/fail)**

- Bitrex, Saccharin.
*Use OSHA approved
fit test protocols
(Appendix A to §
1910.134: Fit Testing
Procedures)*

*Available on OSHA.gov
website*




OSHA Fit Test Guidelines (in part)

Important Points to know prior to testing!

- ▶ The person being tested shall conduct a user seal check (fit check).
- ▶ A valid test can not be conducted if there is any hair growth between the skin and the facepiece sealing surface, such as stubble beard growth, beard, mustache or sideburns which cross the respirator sealing surface. Any type of apparel which interferes with a satisfactory fit shall be altered or removed.
- ▶ ***If the person being tested exhibits difficulty in breathing during the tests, she or he shall be referred to a licensed health care professional to determine whether the test subject can wear a respirator.***
- ▶ The fit test should be performed while the person being tested is wearing any safety equipment that may be worn during actual respirator use.

<https://www.osha.gov/laws-regs/regulations/standardnumber/1910/1910.134AppA>

<https://www.cdc.gov/niosh/npptl/pdfs/FacialHairWmask11282017-508.pdf>



Fit Testing

Before an employee uses any respirator with a tight-fitting face piece, the employee must be fit tested with the same make, model, style, and size of respirator that will be used.



**No respirator fits everyone
and every face -
if it doesn't fit...it doesn't
work!**

Fit Test Solutions for Qualitative Fit Testing

<https://www.osha.gov/laws-regs/regulations/standardnumber/1910/1910.134AppA>

Saccharin Solution – most common

Bitter solution (Bitrex)

- The Bitrex™ (Denatonium benzoate) solution aerosol QLFT protocol uses the same saccharin test protocol because that protocol is widely accepted.
- Note: Bitrex is routinely used as a taste aversion agent in household liquids to help protect children - and is endorsed by the American Medical Association, the National Safety Council, and the American Association of Poison Control Centers.

➤ **Irritant Smoke –**

- Not as common
- Use HEPA (P100) filters
- Only use in well ventilated area
- Not as subjective as saccharin or bitter

➤ **Isoamyl Acetate - Banana Oil**

- Use organic vapor cartridges
- Not as common
- Depends on sense of smell

Translate Medical Evaluation to OSHA Exempt Agriculture

▶ **Best Management Practice for Clinicians**

- ▶ Use OSHA medical evaluation guidelines
- ▶ OSHA Standard 29 CFR 1910.134
- ▶ Ask questions related to work environment, general health, history of respiratory or cardiac conditions
- ▶ Check blood pressure and pulse
- ▶ In addition to basic fit test, observe ability to go around corners, ascend and descend steps, move from high level light to low level light
- ▶ Spend adequate time in care and cleaning education

https://www.cdc.gov/niosh/npptl/topics/respirators/disp_part/respsource3medical.html

Who can perform a Fit Test?





- ▶ **Can anyone perform a respirator fit test?**
- ▶ If the person selected to **perform** the **fit test can** follow the procedure properly and documents the results, it should be considered compliant. It would be wise to make sure the person you designate to perform the fit is **familiar with – and adheres to - the OSHA standards** for compliance.
- ▶ **How to become a Fit Tester:**
 - ▶ Work with someone with respirator fit testing experience
 - ▶ Purchase a fit test kit – step by step instructions are included
 - ▶ Take an online educational course

<https://www.xosafety.academy/courses/performing-respirator-fit-tests>

https://www.3m.com/3M/en_US/safety-centers-of-expertise-us/respiratory-protection/fit-testing/

Pesticide Applicators

- ▶ **Are pesticide applicators and handlers required to wear respirators?**
- ▶ Pesticide applicators and handlers are required to wear respirators if the pesticide label calls for them to do so. Therefore, applicators and handlers need to have had respirator medical evaluation and respirator fit test prior to working with pesticides.
- ▶ **Which pesticides require respirators?**
- ▶ Read and make sure you understand the labels of the pesticides you may use. The label of the pesticides will state if a respirator is required and is based on class of pesticide and application method (drench, spray etc.).
- ▶ <https://fruitgrowersnews.com/article/worker-protection-standard-respirator-rules-updated/>
<http://www.pesticideresources.org/wps/hosted/PERC-WPS-Respirator-Guide.pdf>
<http://pesticideresources.org/med/resources/rmpp.html>



So...what about Availability?

- during a pandemic!



<https://www.cdc.gov/niosh/topics/hcwcontrols/recommendedguidanceextuse.html>

Counterfeit respiratory protection:

<https://www.cdc.gov/niosh/npptl/usernotices/counterfeitResp.html>

What about pesticide respiratory protection availability?

- *during a pandemic!*



- Pesticides may **not** be applied without the label-required PPE.
 - ***There is no substitute for reading the label!***
 - Home-made masks are **not** sufficient substitutes for label-required respirators/masks.
 - **No** exemption or relaxation of the requirements has been made by EPA.
 - Users may need to select alternative products or practices, if the required PPE is **not** available.
- There is a lead educator in each state, focused on Pesticide Safety Education Programs (PSEP). <https://npsec.us/psep-contacts>

What about availability?

- *during a pandemic!*



OSHA Interim Enforcement Response Plan

<https://www.oshalawblog.com/2020/04/articles/osha-interim-enforcement-response-plan-for-covid-19-matters/>


EPA Temporary Guidance on Respiratory Protection in Agriculture

<https://www.epa.gov/pesticides/epa-releases-temporary-guidance-respiratory-protection-agricultural-pesticide-handlers>

Reviewed June 1, 2020

Safety Considerations – What Can We Do?

Engineer dangers out of environment if possible!



Personal protective equipment - may be only feasible solution
The last line of defense!


Know Your Exposures – There are Many!



Gas	Health Effects		
	Low	Medium	High
Hydrogen Sulfide (H₂S)	2-20 ppm: nausea, headache, dizziness	100-300 ppm: altered breathing, fluid in lung	500-700 ppm: collapse, death
Methane (CH₄, LEL)	5-15% (10,000 ppm): asphyxiant	< 1% (10,000 ppm): No known toxicity	< 5% (50,000 ppm): Explosive
Ammonia (NH₃)	5-20 ppm: odor, eye irritation	20-50 ppm: Moderate eye and upper respiratory tract irritation	2500 ppm: chemical pneumonitis, edema, cyanosis, death
Carbon Dioxide (CO₂)	600-1000 ppm: odor, muscle stiffness	1000-2500 ppm: general drowsiness; 5000 ppm: 8-hr maximum	30,000 ppm (3%): increased pulse rate, nausea, impairment
Carbon Monoxide (CO)	< 9 ppm: comfortable living concentration (35 ppm = 8-hr allowable)	200 ppm: headache, dizziness, nausea in 35 ppm = 8-hr 2 hours	400 ppm: life threatening in 3 hours

Children, elderly, pregnant women, and other risk factors: CO, H₂S, and NH₃. The concentrations are advisory only as "no-harm."





When wearing a respirator – the fit is the key

- ▶ All respirators - N95 and above – must fit tight to the face and allow no air leaks
- ▶ If worn correctly, a respirator will filter out 95% or more of aerosol/particulates
- ▶ Respirators must:
 - ▶ be donned correctly
 - ▶ have a fit check performed each time it is worn
 - ▶ be removed (doffed) properly





Image courtesy of the U.S. EPA

Particle Size

Seal Check (Fit check) Procedure: an employee (wearer) responsibility – and is done each time a respirator is used

POSITIVE PRESSURE CHECK



NEGATIVE PRESSURE CHECK



Very Important!

Seven Steps to Correctly Wear a Respirator at Work

Following these simple steps will help you properly put on and take off your respirator, and keep you and everyone else safe.

1 Wash Your Hands




Wash your hands with soap and water or alcohol-based hand rubs containing at least 60% alcohol.

2 Inspect the Respirator



Inspect the respirator for damage. If it appears damaged or damp, do not use it.

3 Put on the Respirator



Clip the respirator in your head with the noseclip at your fingertips and the straps hanging below your head.



Cover your mouth and nose with the respirator and make sure there are no gaps (e.g., facial hair, hair, and glasses) between your face and the respirator.



Place the strap over your head and rest at the top back of your head. If you have a second strap, place the bottom strap around your neck and below your ears. Do not crisscross straps.



If your respirator has a metal nose clip, use your fingertips from both hands to mold the nose area to the shape of your nose.

4 Adjust the Respirator



Place both hands over the respirator. Inhale quickly and then exhale. If you feel leakage from the nose, readjust the noseclip; if leakage from the respirator edges, readjust the straps.



Repeat until you get a proper seal. If you can't get a proper seal, try another respirator.



Avoid touching the respirator while using it. If you do, wash your hands.

Note: If you reuse your respirator, wear gloves when inspecting and putting on the respirator. Avoid touching your face (including your eyes, nose, and mouth) during the process.

6 Remove the Respirator



Wash your hands.



Remove the respirator from behind. Do not touch the front.



If the respirator does not need to be reused because of supply shortages, discard it in a closed-bin waste receptacle. Wash your hands.

For more information, see the quick video, "Putting On and Taking Off a Mask".



Helping You Wear it Right

Wearing Your Filtering Facepiece Respirator

1



Place the respirator over your nose and mouth. Be sure the metal nose clip is on top. With models 8210 or 07048, pre-stretch the straps before wearing.

2



Pull the top strap over your head until it rests on the crown of your head above your ears.

3



Pull the bottom strap over your head until it rests just below your ears.

4



Using both hands starting at the top, mold the metal nose clip around your nose to achieve a secure seal.



Filtering Facepiece (non-valved) Respirator

Filtering Facepiece (non-valved) Respirator

Check the Seal of Your Filtering Facepiece Respirator Each Time You Don the Respirator.



Positive Pressure User Seal Check

For Non-Valved Respirators

Place both hands completely over the respirator and *exhale*. The respirator should bulge slightly. If air leaks between the face and faceseal of the respirator, reposition it and readjust the nose clip for a more secure seal. If you cannot achieve a proper seal, **do not** enter the contaminated area. See your supervisor.



Negative Pressure User Seal Check

For Valved Respirators

Place both hands over the respirator and *inhale* sharply. The respirator should collapse slightly. If air leaks between the face and faceseal of the respirator, reposition it and readjust the nose clip for a more secure seal. If you cannot achieve a proper seal, **do not** enter the contaminated area. See your supervisor.

WARNING

Respiratory Products

These respiratory products protect against certain airborne contaminants. Before use, the wearer must read and understand the User Instructions provided as a part of the product packaging. A written respiratory protection program must be implemented meeting all the requirements of OSHA 1910.134 including training, fit testing and medical evaluation. In Canada, CSA Standard Z94.4 requirements must be met. Respiratory products of the equivalent protection, if appropriate, **failure may result in sickness or death**. For proper use, see packaging instructions, supervisor, or call 3M Chemical Technical Service in USA at 1-800-243-6030 and in Canada at 1-800-387-4434.

General Office:
3M Center
St. Paul, MN 55144-1000

Please recycle. Printed in U.S.A.
© 2014 3M. All rights reserved.
70-0707-0700-0



OSHA Occupational
Safety and Health
Administration
www.osha.gov

1-800-321-OSHA (6742)
TTY 1-877-889-5627

OSHA 4015-02-0000

<https://www.osha.gov/Publications/OSHA4015.pdf>

<https://www.youtube.com/watch?v=oU4stQgCtV8>

Filtering Facepiece Fit Test

Fit Testing is required if:

- ▶ your employer requires the use of respirators.
- ▶ you (the employee) are exposed above the OSHA permissible exposure levels (PELs) of particulate or aerosol.



*If workers wear these respirators
on a voluntary basis NO fit testing
is required.*



Type of Fit Testing

- Quantitative

- Portacount test

- Qualitative (pass/fail)

- Bitrex, Saccharin.
*Use OSHA approved
fit test protocols
(Appendix A to §
1910.134: Fit Testing
Procedures)*

*Available on OSHA.gov
website*




OSHA Fit Test Guidelines (in part)

Important Points to know prior to testing!

- ▶ The person being tested shall conduct a user seal check (fit check).
- ▶ A valid test can not be conducted if there is any hair growth between the skin and the facepiece sealing surface, such as stubble beard growth, beard, mustache or sideburns which cross the respirator sealing surface. Any type of apparel which interferes with a satisfactory fit shall be altered or removed.
- ▶ ***If the person being tested exhibits difficulty in breathing during the tests, she or he shall be referred to a licensed health care professional to determine whether the test subject can wear a respirator.***
- ▶ The fit test should be performed while the person being tested is wearing any safety equipment that may be worn during actual respirator use.

<https://www.osha.gov/laws-regs/regulations/standardnumber/1910/1910.134AppA>

<https://www.cdc.gov/niosh/npptl/pdfs/FacialHairWmask11282017-508.pdf>



Fit Testing

Before an employee uses any respirator with a tight-fitting face piece, the employee must be fit tested with the same make, model, style, and size of respirator that will be used.



**No respirator fits everyone
and every face -
if it doesn't fit...it doesn't
work!**

Fit Test Solutions for Qualitative Fit Testing

<https://www.osha.gov/laws-regs/regulations/standardnumber/1910/1910.134AppA>

Saccharin Solution – most common

Bitter solution (Bitrex)

- The Bitrex™ (Denatonium benzoate) solution aerosol QLFT protocol uses the same saccharin test protocol because that protocol is widely accepted.
- Note: Bitrex is routinely used as a taste aversion agent in household liquids to help protect children - and is endorsed by the American Medical Association, the National Safety Council, and the American Association of Poison Control Centers.

➤ **Irritant Smoke –**

- Not as common
- Use HEPA (P100) filters
- Only use in well ventilated area
- Not as subjective as saccharin or bitter

➤ **Isoamyl Acetate - Banana Oil**

- Use organic vapor cartridges
- Not as common
- Depends on sense of smell

Translate Medical Evaluation to OSHA Exempt Agriculture

▶ **Best Management Practice for Clinicians**

- ▶ Use OSHA medical evaluation guidelines
- ▶ OSHA Standard 29 CFR 1910.134
- ▶ Ask questions related to work environment, general health, history of respiratory or cardiac conditions
- ▶ Check blood pressure and pulse
- ▶ In addition to basic fit test, observe ability to go around corners, ascend and descend steps, move from high level light to low level light
- ▶ Spend adequate time in care and cleaning education

https://www.cdc.gov/niosh/npptl/topics/respirators/disp_part/respsource3medical.html

Who can perform a Fit Test?





- ▶ **Can anyone perform a respirator fit test?**
- ▶ If the person selected to **perform** the **fit test can** follow the procedure properly and documents the results, it should be considered compliant. It would be wise to make sure the person you designate to perform the fit is **familiar with – and adheres to - the OSHA standards** for compliance.
- ▶ **How to become a Fit Tester:**
 - ▶ Work with someone with respirator fit testing experience
 - ▶ Purchase a fit test kit – step by step instructions are included
 - ▶ Take an online educational course

<https://www.xosafety.academy/courses/performing-respirator-fit-tests>

https://www.3m.com/3M/en_US/safety-centers-of-expertise-us/respiratory-protection/fit-testing/

Pesticide Applicators

- ▶ **Are pesticide applicators and handlers required to wear respirators?**
- ▶ Pesticide applicators and handlers are required to wear respirators if the pesticide label calls for them to do so. Therefore, applicators and handlers need to have had respirator medical evaluation and respirator fit test prior to working with pesticides.
- ▶ **Which pesticides require respirators?**
- ▶ Read and make sure you understand the labels of the pesticides you may use. The label of the pesticides will state if a respirator is required and is based on class of pesticide and application method (drench, spray etc.).
- ▶ <https://fruitgrowersnews.com/article/worker-protection-standard-respirator-rules-updated/>
<http://www.pesticideresources.org/wps/hosted/PERC-WPS-Respirator-Guide.pdf>
<http://pesticideresources.org/med/resources/rmppp.html>



So...what about Availability?

- during a pandemic!



<https://www.cdc.gov/niosh/topics/hcwcontrols/recommendedguidanceextuse.html>

Counterfeit respiratory protection:

<https://www.cdc.gov/niosh/npptl/usernotices/counterfeitResp.html>

What about pesticide respiratory protection availability?

- *during a pandemic!*



- Pesticides may **not** be applied without the label-required PPE.
 - ***There is no substitute for reading the label!***
 - Home-made masks are **not** sufficient substitutes for label-required respirators/masks.
 - **No** exemption or relaxation of the requirements has been made by EPA.
 - Users may need to select alternative products or practices, if the required PPE is **not** available.
- There is a lead educator in each state, focused on Pesticide Safety Education Programs (PSEP). <https://npsec.us/psep-contacts>

What about availability?

- *during a pandemic!*



OSHA Interim Enforcement Response Plan

<https://www.oshalawblog.com/2020/04/articles/osha-interim-enforcement-response-plan-for-covid-19-matters/>

EPA Temporary Guidance on Respiratory Protection in Agriculture

<https://www.epa.gov/pesticides/epa-releases-temporary-guidance-respiratory-protection-agricultural-pesticide-handlers>

Reviewed June 1, 2020